

Targeted Isolation of Tumor-Derived Extracellular Vesicles Using a Cocktail of Magnetic Nanospheres for Minimally Invasive Pancreatic Cancer Diagnosis

1. INTRODUCTION AND BACKGROUND

Pancreatic Cancer – The King of Cancers


Epidemiology	Prognostics	Diagnosis and treatment
4 th leading cause of cancer-related deaths in U.S. today	1 st mortality rate	87% diagnosed at advanced stages
2 nd leading cause of cancer-related deaths in U.S. by 2030	4 months median overall survival	< 20% eligible for surgery
10 th in incidence in U.S. today		

Tables created by finalist in PowerPoint using data from Stoop, T. F., et. al. Pancreatic cancer. *The Lancet* 405, 1182-1202, (2025).

Current Challenges in Pancreatic Cancer Diagnosis

- **Medical imaging** is widely used for cancer detection but is limited by *high cost, low sensitivity, radiation exposure risks, and infrequent applicability*
 - **Tissue biopsy** enables molecular diagnosis, yet it is *invasive, expensive, prone to sampling bias, and often inaccessible for pancreatic cancer*
 - **Circulating free DNA-based liquid biopsy** offers a minimally invasive alternative but *lacks sensitivity for early cancer diagnosis and has limited multiplexing capability for RNA and protein analysis*
- ➔ To win the battle against pancreatic cancer, it is urgent to develop innovative and effective diagnostic methods for early diagnosis, treatment monitoring, and recurrence surveillance.

2. APPROACH: EV ISOLATION + ddPCR



3. RESULTS

Comparison with Standard Method (Exploratory Cohort)

KRAS mutations were detected in all 20 pancreatic cancer patients, while 9 of 10 healthy donors were mutation-negative.


Sensitivity and specificity

	Patients	Healthy donors
Test positive	TP = 20	FP = 1
Test negative	FN = 0	TN = 9

Sensitivity = 100%

Specificity = 90%


Discrimination



EV-DNA amount

Antibody cocktail-based EV isolation yielded significantly less EV DNA than total EV isolation (Wilcoxon matched-pairs signed rank test $p = 0.0459$).


Figures created by finalist using Prism, 2025



Clinical sensitivity

The antibody cocktail-based EV isolation achieved **100%** KRAS detection versus **73.9%** with total EV isolation, with significantly higher detectability (Wilcoxon matched-pairs signed rank test $p = 0.0422$).

Figures created by finalist using Prism, 2025



Impact of Cancer Stage on EV-based Clinical Diagnosis

Pooling training and validation cohorts, PDAC patients were grouped as local/early (stages I–II) or advanced (stages III–IV) stage disease.

EV-DNA amount




Figure created by finalist using Prism, 2025

EV-DNA yields were higher in advanced versus local disease ($p = 0.0444$, Mann-Whitney).

Early-stage sensitivity, specificity and discrimination

Sensitivity and specificity

	Patients	Healthy donors
Test positive	TP = 15	FP = 2
Test negative	FN = 4	TN = 16

Sensitivity = 78.9%

Specificity = 88.9%

Table created by finalist in PowerPoint, 2025

Discrimination




Figure created by finalist using Prism, 2025

KRAS mutation MAF




Figure created by finalist using Prism, 2025

KRAS MAF was higher in advanced versus local disease ($p = 0.0469$, Mann-Whitney).

Late-stage sensitivity, specificity and discrimination

Sensitivity and specificity

	Patients	Healthy donors
Test positive	TP = 25	FP = 2
Test negative	FN = 2	TN = 16

Sensitivity = 92.6%

Specificity = 88.9%

Table created by finalist in PowerPoint, 2025

Discrimination




Figure created by finalist using Prism, 2025

Evaluation in Validation Cohort

KRAS mutations were detected in 20 pancreatic cancer patients (6 false negatives), while 7 of 8 healthy donors were mutation-negative.

Sensitivity and specificity

	Patients	Healthy donors
Test positive	TP = 20	FP = 1
Test negative	FN = 6	TN = 7

Sensitivity = 76.9%

Specificity = 87.5%

Table created by finalist in PowerPoint, 2025

Discrimination




Figure created by finalist using Prism, 2025

Overall sensitivity, specificity and discrimination

Sensitivity and specificity

	Patients	Healthy donors
Test positive	TP = 40	FP = 2
Test negative	FN = 6	TN = 16

Sensitivity = 87.0%

Specificity = 88.9%

Table created by finalist in PowerPoint, 2025

Discrimination




Figure created by finalist using Prism, 2025

4. CONCLUSIONS

- ✓ **87.0% sensitivity and 88.9% specificity** for 46 pancreatic cancer patients and 18 healthy donors.
- ✓ **Compared to total EV isolation**, cancer-specific EV increased clinical sensitivity (100% vs 73.9%).
- ✓ **Overall method** offers a highly sensitive, minimally invasive approach for molecular diagnosis.

KEY REFERENCES

1. Kalluri, R., and LeBleu, V. S. (2020), *Science*, 367(6478), eaau6977.
2. Shao, H., et. al. *Chemical Reviews*, 118(4), 1917-1950.
3. Weinberg, R. A. (2013). *Biology of Cancer*, Second Edition. Garland Science.